



Mindray Hematology Solution for **Oncology Hospitals**

Advancing Hematology Solutions to Elevate Diagnosis
and Treatment of Cancer Care



Lab's Challenges

- Detection capabilities for immature and abnormal cells.
- Handling challenging samples such as Low-PLT or abnormal cells in peripheral blood for prognosis.
- Lack of pathologists during night shifts.



Lab's Needs

- A reliable analyzer for detecting immature and abnormal cells by cross-referencing morphological results with hematology data in the presence of abnormal samples.
- Solutions for samples with low WBC and PLT counts.
- Remote review of critical samples anytime, anywhere.

Mindray BC-700 Series / CAL 8000 with MC-80 for Medium and Large Samples Laboratories

Hematology-Oncology



BC-700 Series




CAL 8000 with Digital Cell Morphology System MC-80

With Mindray 3D SF Cube Technology for Abnormal Cells Detection

DIFF Scattergram




- ❶ Abnormal lymphocytes
- ❷ Immature granulocytes
- ❸ Monocytes
- ❹ Lymphocytes




- ❺ Neutrophils and basophils
- ❻ Eosinophils
- ❼ Blood ghost

WNB Scattergram



- ❶ Nucleated red blood cells
- ❷ Basophils
- ❸ White blood cells
- ❹ Blood ghost

RET Scattergram



- ❶ Mature red blood cells
- ❷ Low-fluorescence reticulocytes
- ❸ Medium-fluorescence reticulocytes
- ❹ High-fluorescence reticulocytes
- ❺ Optical platelets

Immature and Abnormal Cells Detection for Early Diagnosis and Prognosis

NRBC, Blast and Immature Granulocytes (IMGs) Detection

The presence of NRBCs, blast cells, and IMGs in peripheral blood samples can serve as important diagnostic indicators of certain types of leukemia, lymphoma, and other hematological malignancies.

NRBC


Reportable parameter by direct measurement for accurate WBC+DIFF results

Blast cells

High sensitivity to avoid missed diagnosis

IMGs

Actual IMG count and flag to reduce review rate



Mindray's Performance in Immature and Abnormal Cells Detection

	All flags	All blast	IG	NRBC
TP	359	122	240	133
TN	198	483	361	489
FP	68	32	31	14
FN	15	3	8	4
Sensitivity	96.0	97.6	96.8	97.1
Specificity	74.4	93.8	92.1	97.2
PPV	84.1	79.2	88.6	90.5
NPV	93.0	99.4	97.8	99.2
Efficiency	87.0	94.5	93.9	97.2

Abbreviations: FN, false negatives; FP, false positives; IG, immature granulocytes; NPV, predictive value of negatives; NRBC, nucleated red blood cells; PPV, predictive value of positives; TN, true negatives; TP, true positives.

Zini, G., Cantelli, F., Scavone, F., Barbagallo, O., & Ciminello, A. (2020). Hematological performance of a last generation automated blood cell counter: The Mindray BC-6800 Plus. *International Journal of Laboratory Hematology*, 42(4), 439–449. <https://doi.org/10.1111/ijlh.13218>

Morphological Blood Cells Confirmation for Early Diagnosis

Cross Checking with Hematology Analyzer for More Accurate Result



Mindray's labXpert system integrates morphology results with hematology data, enabling precise diagnostic assessments in healthcare settings.

Advanced imaging technology, unified data management, and streamlined reporting facilitate efficient workflow, enhance quality control, and support evidence-based decision-making. This integration optimizes patient care and diagnostic accuracy in oncology and hematology practices.

- 1 Early Detection
- 4 Prognostic Assessment
- 2 Disease Monitoring
- 5 Treatment Guidance
- 3 Subtyping and Classification

Research on White Blood Cell (WBC)

Compared with the expert-reviewed results of 146,538 cells, the AI pre-classification results have an accuracy rate as high as 97.8%, fully demonstrating that AI has extremely high pre-classification capabilities.




AI Pre-Classification

VS




Morphologist Reviewed Result
(Morphologist review was conducted on the AI pre-classification results)




Blast Cells

Accuracy: **89.58%**
Sensitivity: **98.38%**
Specificity: **99.40%**




Abnormal Lymphocytes

Accuracy: **68.44%**
Sensitivity: **97.26%**
Specificity: **99.63%**



Abnormal Promyelocytes


Accuracy: **89.96%**
Sensitivity: **80.00%**
Specificity: **99.92%**




Plasma Cells

Accuracy: **96.61%**
Sensitivity: **73.08%**
Specificity: **100.00%**


Blood Cell Images Captured from MC-80



NRBC



Blast Cells



IMG

Handling Low Count WBC & PLT Samples Without Human Intervention

Automatic Multiple Counting of Low WBC & PLT

Low WBC

- Inaccurate WBC and DIFF results
- Missed diagnosis of blast cells leading to delay diagnosis & treatment




- 1 Accurate WBC & DIFF results are needed
- 2 Confirm the **blast cells** for early leukemia diagnosis
- 3 Prompt the next step for treatment

3x

Number of Particles

Ensure accurate determination of low-value cells



Low PLT

- Inaccurate PLT result
- Pseudo-thrombocytopenia leading to delay surgery & wrong medication




- 1 Accurate PLT result is needed
- 2 Indicated **giant PLT/PLT aggregation**
- 3 Free from **interferences** e.g. RBC fragments and small RBC

8x

Number of Particles

Ensure accurate determination of low-value cells



Multiple Counting for Low-WBC Samples



For leucopenia samples ($WBC < 1.5 \times 10^9/L$) or with capillary blood samples, we automatically extend the count time to increase the number of cells counted. This process enhances the accuracy and reliability of WBC results.

PLT-O with Multiple Counting for Low-Count Platelet Samples

To begin with, PLT-I value obtained from impedance channel is firstly compared with a default value ($50 \times 10^9/L$). If this is lower than the cutoff value, the analyzer can automatically prolong counting time up to 8-fold in order to collect more PLT particles for further analysis. In addition, Auto 8xPLT-O Counting Technology can eliminate other interference factors (e.g. fragmented RBC/WBC) which are often easily miscounted as PLT by the impedance channel. No additional sampling, no manual intervention, no additional channels and reagents required. 8xPLT-O counting is both efficient and effective in counting low PLT accurately.




Total PLT Solution for Better Cancer Patients Monitoring for Treatment

Accurate PLT Count with Images Captured from MC-80



Total PLT Solution for Anti-Interference and Accurate PLT Results



Chemotherapy is a cornerstone in cancer treatment, but its side effects can pose challenges for patients, the drugs can directly affect the bone marrow, impairing platelet production.


Clinical Case

A patient with systematic lupus erythematosus revisited the hospital and got checked by Mindray BC-700 Series hematology analyzer.

Cbc Results from BC-700 Series


MCV	98.8 fL
MPV	17.0 fL
P-LCR	75.2%
PLT-I	52x 10 ⁹ /L
PLT-H	89x 10 ⁹ /L
PLT-O	91x 10 ⁹ /L

Microscopic examination revealed evident large PLTs (arrows) having a similar size to the RBCs in multiple high-power lens.



Method	PLT-I	PLT-H	PLT-O	PLT-M	Flow Cytometry
PLT(x 10 ⁹ /L)	52	89	91	92	90


The PLT count measured by flow cytometry was 90x10⁹/L, which was consistent with that by PLT-H. In this case, the MPV and P-LCR values increased significantly. The tail of the PLT histogram was elevated, suggesting that large PLTs may be present in this sample, which would cause a pseudo-thrombocytopenia. The information obtained from the DIFF channel can ensure the accurate detection of large PLTs, avoiding unnecessary blood transfusion.



Revolutionize Lab Workflow with labXpert System

Real-Time Integration for Remote Diagnostics and Enhanced Clinical Expertise Utilization

Elevate your laboratory capabilities with labXpert software, enabling real-time integration and sharing of testing results. Experience the benefits of remote diagnostics and optimized utilization of clinical professional expertise, a valuable resource in high demand. Streamline operations, enhance efficiency, and drive better outcomes in your lab with labXpert system.



Affiliated Hospital A

Affiliated Hospital B

Remote Consultation


Head office: FusionGreen Healthcare Pvt Ltd
No 18/1, Portion A, Ground Floor, 53rd Street, Ashok Nagar,
Chennai-600 083, Tamil Nadu.

Branch Offices: Coimbatore, Trichy, Madurai & Vellore



www.fusiongreen.in

Marketed By :



+ 91 96771 44414, 99401 41956